

無水エタノール希釈割合（秋田放射能測定室「べぐれでねが」作成）

・エタノール濃度70%

作成する消毒液(ml)	無水エタノール(ml)	精製水(ml)
10	7	3
20	14	6
30	21	9
40	28	12
50	35	15
60	42	18
70	49	21
80	56	24
90	63	27
100	70	30
110	77	33
120	84	36
130	91	39
140	98	42
150	105	45
160	112	48
170	119	51
180	126	54
190	133	57
200	140	60
210	147	63
220	154	66
230	161	69
240	168	72
250	175	75
260	182	78
270	189	81
280	196	84
290	203	87
300	210	90
310	217	93
320	224	96
330	231	99
340	238	102
350	245	105
360	252	108
370	259	111
380	266	114
390	273	117
400	280	120
410	287	123
420	294	126
430	301	129
440	308	132
450	315	135
460	322	138
470	329	141
480	336	144
490	343	147
500	350	150
510	357	153
520	364	156
530	371	159
540	378	162
550	385	165
560	392	168
570	399	171
580	406	174
590	413	177
600	420	180
610	427	183
620	434	186
630	441	189
640	448	192
650	455	195
660	462	198
670	469	201
680	476	204
690	483	207
700	490	210
710	497	213
720	504	216

・エタノール濃度80%

作成する消毒液(ml)	無水エタノール(ml)	精製水(ml)
10	8	2
20	16	4
30	24	6
40	32	8
50	40	10
60	48	12
70	56	14
80	64	16
90	72	18
100	80	20
110	88	22
120	96	24
130	104	26
140	112	28
150	120	30
160	128	32
170	136	34
180	144	36
190	152	38
200	160	40
210	168	42
220	176	44
230	184	46
240	192	48
250	200	50
260	208	52
270	216	54
280	224	56
290	232	58
300	240	60
310	248	62
320	256	64
330	264	66
340	272	68
350	280	70
360	288	72
370	296	74
380	304	76
390	312	78
400	320	80
410	328	82
420	336	84
430	344	86
440	352	88
450	360	90
460	368	92
470	376	94
480	384	96
490	392	98
500	400	100
510	408	102
520	416	104
530	424	106
540	432	108
550	440	110
560	448	112
570	456	114
580	464	116
590	472	118
600	480	120
610	488	122
620	496	124
630	504	126
640	512	128
650	520	130
660	528	132
670	536	134
680	544	136
690	552	138
700	560	140
710	568	142
720	576	144